

GENESY 3000

ELECTROTHERAPY 4 CHANNELS | 423 PROGRAMS

TOP PROFESSIONAL FOUR-CHANNEL PRODUCT, USED IN THE MEDICAL FIELD FOR THE TREATMENT OF THE MOST COMMON PATHOLOGIES OF THE MUSCLE-SKELETICAL SYSTEM.

236 ELECTROTHERAPY MENU (EMS/RUSSIAN) | 53 SPORT | 33 PAIN | 29 FITNESS | 23 MICROCURRENTS | 16 BEAUTY | 11 REHAB | 9 DENERVATED | 9 INTERFERENTIALS | 3 INCONTINENCE | 1 IONOPHORESIS


TECHNICAL FEATURES

Display
Visible area size 2,6"
Channels
4 independent (8 electrodes)
Frequency
0,3-150Hz *
Pulse amplitude
40-450µs *
Power
0-120mA per channel
Power supply
- Mains - Rechargeable batteries
Size
mm 170x220x60
Weight
gr 1152


* according to the type of current

EQUIPMENT

- 1 Bag
- 1 GENESY 3000 stimulator
- 4 Cables for electrode connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 2 Iono electrodes + 2 elastic bands
- 1 Power supply unit
- 1 Operating manual



CE
0476



GLOBUS
ITALIAN EXCELLENCE

GENESY 3000

ELECTROTHERAPY 4 CHANNELS | 423 PROGRAMS

PROGRAMS LIST

REHAB

Swollen ankles

Atrophy recovery

Hemiplegy-upper limbs

Hemiplegy-lower limbs

Recovery after ACL surgery

Functional recovery

Ankle re-education

Leg re-education

Muscle Spasms

Shoulder subluxation prevention

Vastus medialis reinforcement

Motor point pen

Agonist/Antagonist

Muscle reinforcement

PAIN

High-frequency antalgic tens

Low frequency antalgic tens

Endorphinic tens

Knee pain

Post-surgical pain

Chronic pain

Shoulder pain

Carpal tunnel

Trapezius pain

Osteoarthritis

Knee osteoarthritis

Rotator cuff tendinitis

Muscle pain

Nerve compression

Spinal compression pain

Spinal osteoarthritis

Spinal osteoporosis

Hip osteoarthritis

Ankle osteoarthritis

Muscle tendon injury pain

Post-surgical pain

Sciatica

Chronic lumbago

Cervical pain

Acute pain inguinal hernia

Bursitis-tendinitis

Muscle injuries

Bone fractures

Whiplash

Epicondylitis

Menstrual pain

Conventional antalgic tens

Modulated antalgic tens

IONOPHORESIS

Ionophoresis

INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

MICROCURRENTS

Epicondylitis

Scapulothoracic periarthritis

Muscle restoration

Contusion

Edema

Skin ulcer

Sciatica

Lumbago

Brachial neuralgia

Acute pain

Articular pain

Stiff neck

Whiplash

Cervical spondylosis

Shoulder sprain

Carpal tunnel

Knee sprain

Knee osteoarthritis

Ankle sprain

Achilles tendon inflammation

Patella tendon inflammation

Rotator cuff inflammation

Tendon inflammation

DENERVATED

Triangular 1

Triangular 2

Triangular 3

Trapezoidal 1

Trapezoidal 2

Trapezoidal 3

Rectangular 1

Rectangular 2

Rectangular 3

INTERFERENTIALS

Knee osteoarthritis

Acute lumbago

Cervical pain

Chronic lumbago

Frozen shoulder

Knee post-surgery pain

Pain vertebral fractures

Knee post-surgery edema

Fibromyalgia

SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Aerobic endurance

Reactivity

Post-competition recovery

Decontracting

Hypertrophy

FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

Aerobic fitness

Cramp prevention

BEAUTY

Drainage

Lipolysis

Toning massage

Skin tone improvement

Connective massage

Post-pregnancy drainage

Post-pregnancy lipolysis

Post-pregnancy firming

ELECTROTHERAPY

The Electrotherapy program list includes 192 EMS and 44 Russian(KOTZ) parameters combinations