

# Circle Line



## Perfection – from professionals for professionals

The concentric-excentric operating dual-function devices for two muscle groups each stand out thanks to their simple engineering and quality workmanship. Weights to increase resistance were not used. Instead, CircleLine works using an own made hydraulic cylinder system in which the resistance is adjusted using an easy-to-use scaled dial. Depending on the device, you can choose whether the level of resistance should be adjusted together or variably for the two directions of movement. This guarantees a smoother and more uniform motion resistance, reduces the risk of injury and overloading, and enables an effective workout in a very short time, since the concentric-excentric system acts on two muscle groups simultaneously.



■ See inside back page for upholstery colour chart



### ABDOMINAL/BACK TRAINER

Combination machine with 2-way cylinder

Item No. HZ 4  
Differentiable resistance for agonist & antagonist

Dimensions:  
L x W x H / 79 cm x 69 cm x 127 cm  
Weight: 44 kg

Option  
Item No. 10009536  
Designed as a medical device - MPG compliant  
Upholstery quality - Palma



### LEG CURL/EXTENSION MACHINE

Combination machine with 2-way cylinder

Item No. HZ 3  
Differentiable resistance for agonist & antagonist

Dimensions:  
L x W x H / 152 cm x 69 cm x 114 cm  
Weight: 63 kg

Option  
Item No. 10009536  
Medical device version - MPG compliant  
Item No. HZ-P  
Upholstery quality - Palma



### UPPER BACK / PECS TRAINER

Combination machine with 2-way cylinder

Item No. HZ 5  
Differentiable resistance for agonist & antagonist

Dimensions:  
L x B x H / 79 cm x 80 cm x 108 cm  
Weight: 49 kg

Option  
Item No. 10009536  
Medical device version - MPG compliant  
Item No. HZ-P  
Upholstery quality - Palma



### SHOULDER LIFT / PULL DOWN

Combination machine with 2-way cylinder

Item No. HZ 9  
Differentiable resistance for agonist & antagonist

Dimensions:  
L x W x H / 92 cm x 69 cm x 126 cm  
Weight: 52 kg

Option  
Item No. 10009536  
Designed as a medical device - MPG compliant  
Item No. HZ-P  
Upholstery quality - Palma



### OUTER/INNER THIGH TRAINER

Combination machine with 2-way cylinder

Item No. HZ 6  
Differentiable resistance for agonist & antagonist

Dimensions:  
L x W x H / 135 cm x 69 cm x 107 cm  
Total weight: 66 kg

Option  
Item No. 10009536  
Designed as a medical device - MPG compliant  
Item No. HZ-P  
Upholstery quality - Palma



### BUTTERFLY/CORE STABILISER

Combination machine with 2-way cylinder

Item No. HZ 8  
Differentiable resistance for agonist & antagonist

Dimensions:  
L x W x H / 104 cm x 69 cm x 113 cm  
Total weight: 67 kg

Option  
Item No. 10009536  
Designed as a medical device - MPG compliant  
Item No. HZ-P  
Upholstery quality - Palma



### LEG PRESS

Item No. HZ 1  
Leg press

Dimensions:  
L x W x H / 154 cm x 69 cm x 123 cm  
Total weight: 71 kg

Option  
Item No. 10009536  
Designed as a medical device - MPG compliant  
Item No. HZ-P  
Upholstery quality - Palma



### CIRCLE COLUMN

With its striking slender design, this low-height column provides CIRCLE Line timekeepers with a valuable accessory for group training.

The LED equipment based on the traffic light principle can be configured to suit your requirements. PATENTED DESIGN

Item No. 10009533  
• Circle column with LED lighting

Dimensions:  
B x D x H / 40 cm x 7 cm x 127 cm